

BioProTT™ FlowMCP Platform (up to 4 channels)



- Transit time ultrasonic flow meter
- Simultaneous measurement of multiple flow channels
- Switch between digital and analog output depending on variant
- Compact space saving design with DIN Rail mounting

Technical Specification

BioProTT™ Variants	FlowMCP 1	FlowMCP 1-a	FlowMCP 2	FlowMCP 2-a	FlowMCP 4
Flow Channels	1	1	2	2	4
Weight	200 g	200 g	240 g	240 g	320 g
Size (HxWxD)	99 x 45 x 115 mm	99 x 45 x 115 mm	99 x 45 x 115 mm	99 x 68 x 115 mm	99 x 90 x 115 mm
Digital Interface	Modbus TCP (RJ-45)	Modbus TCP (RJ-45)	Modbus TCP (RJ-45)	Modbus TCP (RJ-45)	Modbus TCP (RJ-45)
Analog Interface		(0) 4-20mA		(0) 4-20mA	
Housing material		Polyamide			
IP-Code		IP20			
Supply voltage		24 V DC (± 10 %) max 2 A via external supply			
Power consumption		Typically 3 Watt, max. 8 Watt - max. 9 Watt at MCP4			
DIN Rail mounting type		Standard TS35 DIN Rail according to EN 60715, 35 mm x 7.5 mm			
Flow sensor extension cables		1.1 m, IP 50 at panel; w/ D-Sub (flow meter) and round (sensor) connector to allow sensor connection from outside control cabinet			
On-site calibration adjustment		With user adaptable calibration factor (per flow channel)			

Compatible Transducer

Full range of BioProTT™ Clamp-On Transducers	For maximum flow measurement range, accuracy and resolution, please see BioProTT™ Clamp-On Transducer datasheet
--	---

Ambient Conditions

Air pressure	70 to 106 kPa
Operating temperature range	10 to 40 °C (50 to 104 °F)
Storage temperature range	-20 to 45 °C (-4 to 113 °F)
Transport temperature range	-20 to 55 °C (-4 to 131 °F)
Humidity storage, transport and operation	10 to 96 % (non-condensing)

D214-705 BioProTT™ FlowMCP 1-4 - Technical Data Sheet - V6.0 | EMT-10001-T-02-A4

Copyright 2021 PSG®, a Dover company

Authorized PSG® Partner:

Where Innovation Flows



em-tec
BIOPROCESSING

em-tec GmbH
Lerchenberg 20

86923 Finning, Germany

P: +49 8806 9236 0

F: +49 8806 9236 50

E: em-tec-info@psgdover.com

em-tec.de